**REF**

GD5428 00

GLICOTEST HbA1c Control/Controllo

Controls for Glycated Hemoglobin HbA1c assay

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**IVD****INTENDED USE**

This product is for the purpose of the control of the results in the quantitative determination of human hemoglobin A1c (HbA1c) in blood by automated immunoassay.

COMPOSITION

1. **Normal Control – Level 1 (lyophilized)** 1 x 0.5 ml
2. **High Control – Level 2 (lyophilized)** 1 x 0.5 ml

Stabilized hemolysate prepared from human erythrocytes.

Reference values

(%) HbA1c assigned values, and relative acceptable limits, are reported in the attached value sheet and are lot specific. The limits of the test could be used as a guide to detect the accuracy of the of the test procedure.

AUXILIARY EQUIPMENT

- Glicotest HbA1c REF GD5424 00
- 0.5 ml pipette
- Deionized water

PREPARATION

Reconstitute each control using 0.5 ml of deionized water. Gently mix for 10 minutes, or until all material has dissolved.

Subdivide the reconstituted control into aliquots of 0.1 ml and seal tightly.

STORAGE AND STABILITY

Store at 2- 8°C. The controls are stable until expiration date if stored at the indicated conditions, protected from light and heat.

After reconstitution:

The reconstituted controls are stable 30 days at 2-8 °C and 3 months at -20 °C. Avoid freeze-thawing cycle.

PROCEDURE

The lyophilized HbA1c control should be treated in the same manner as specimens, hemolysate procedure included. Follow the instructions for use accompanying the kits and the applications for the instruments.

Test limitation

Improper pipetting, inadequate mixing and erroneous instruments calibration may cause inaccurate results.

PRECAUTIONS IN USE

All reagents contain inactive components such as preservatives (Sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/ECC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.

The use of the laboratory reagents according to good laboratory practice is recommended.

The reagents from human donors have given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle with caution.

Waste Management

Please refer to local legal requirements.

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